

# How to Check Hearing Aid Batteries

Hearing aids batteries need to be replaced regularly for the hearing aids to work.

You can use a battery tester to check how your child's battery is working. There are several types, but they all work in essentially the same way. They have an indicator bar or line that shows you the strength of the battery.

Place the battery on the tester. If a battery is good, the indicator bar or line will show you the strength of the battery. If the battery is not working or is losing strength, it is time to change the battery.

Some hearing aids have an indicator light that flashes to indicate when the battery is charged or needs to be replaced. Ask your audiologist how the light on your child's hearing aid works.