

How to Clean the Earmold

Wax can get into the earmold opening and block the sound. First, remove any wax build up from the earmold opening with the wire loop that came with the hearing aid care kit. You can do this every day, and if needed, you can wash the earmold.

The first step is to detach the earmold from the hearing aid.

Second, soak the earmold in warm, soapy water, liquid dish soap works well. Rinse it and then wipe it dry with a towel. It works well to do this at night so you can let it air dry for several hours.

The third step is to use the air blower to remove moisture from the earmold tubing before re-attaching the earmold to the hearing aid.